

..... **outstanding, novel**

**Miniaturization of HPLC**, such as  
Capillary HPLC (on-line LC/MS) Microbore HPLC requires  
thin columns in order to achieve maximum of sensitivity


..... therefore urgently needed:

Chemically resistant glass microtiterplates for sample storage,  
robotic sample application, fraction collection, combinatorial  
chemistry and many other applications prior to or following  
HPLC.

Deep Well Microtiter **GROM-** Plates  
made from glass


**96-well**

*Deep Well Microtiter Plate*



**384-well**

*Deep Well Microtiter Plate*



**1536-well**

*Deep Well Microtiter Plate*

For more detailed informations see pages 91 and 92

# HPLC technologies:

## Combinatorial Chemistry

requires

short columns to dramatically increase sample throughput (HTS)

for highly efficient and rapid chromatography:

### **NovoGROM** low-dispersion, High Speed HPLC columns



50 x 0,3 mm i.d. capillary High Speed column



50 x 1.0 mm i.d. microbore High Speed column




50 x 4.0 mm i.d. analytical High Speed column



50 x 8,0 mm i.d. semipreparative High Speed column



50 x 20 mm i.d. preparative High Speed column



50 x 50 mm i.d. preparative High Speed column


For more detailed informations see pages 87-90

## Ultra-stable **Zirconia-based** stationary phases outclassing polymeric and other HPLC columns

**GROM Zirconia** phase-packed HPLC columns represent superior alternative to polymeric columns. They possess all of the advantages of polymeric columns with none of the drawbacks and show extreme chemical and thermal stability.


**GROM Zirconia** phase-packed columns show superior efficiency, allow highly reproducible separations particularly for gradient elution, and exhibit selectivity similar to bonded phases for easier method development.

### Analysis of Benzoates by Novel Zirconia Phases



- 1) Methyl benzoate
- 2) Ethyl benzoate
- 3) Propyl benzoate
- 4) Butyl benzoate
- 5) Pentyl benzoate


### Influence of linear flow on HETP



**Column phase:** GROM-Zirconia RPP, 3 µm  
**Column size:** 50 x 2 mm  
**Eluent:** H<sub>2</sub>O / ACN = 65 / 35  
**Flow rate:** 1.0 ml/min  
**Pressure:** 16 MPa  
**Temperature:** RT  
**Detection (UV):** 254 nm (1.2 µl flow cell)  
**Injection:** 2 µl (50-150 µg/ml of each)

## .... outstanding features

- Extraordinary pH (1-14) and unique temperature stability ( $\geq 200^\circ\text{C}$ )
- Superior efficiency
- Exceptional solvent resistance and selectivity very similar to bonded phases
- ➔ At least equivalent to all polymer brands
- ➔ Nearly twice as efficient as most polymer brands
- ➔ Enables easy method developments



### Analysis of


- 1) Resorcinol, 2) Benzonitrile,
- 3) Methyl benzoate, 4) Anisole, 5) Benzene

**Column phase:** Zirconia RP P, 3 µm  
**Column size:** 150 x 4.6 mm  
**Eluent:** ACN / H<sub>2</sub>O = 15 / 85  
**Flow rate:** 1.0 ml/min  
**Pressure:** 18 MPa  
**Temperature:** 50°C  
**Detection (UV):** 254 nm  
**Injection:** 5 µl

# HPLC columns packed with


## Zirconia phases guaranteeing unmatched chemical and thermal stability

### Gradient Separation of n-Alkylbenzenes by Novel Zirconia Phases




**Column phase:** GROM-SIL Zirconia RP P, 3 µm  
**Column size:** 50 x 2 mm  
**Eluent A:** H<sub>2</sub>O  
**B:** ACN  
**Gradient:** 30 - 80% B (0 - 2.5 min)  
**Flow rate:** 0.6 ml/min  
**Pressure:** 10 MPa  
**Temperature:** RT  
**Detection (UV):** 254 nm (1.2 µl flow cell)  
**Injection:** 2 µl (3 mg/ml of each)

### Thermal Stability at 195°C



### Gradient Separation of Benzoates by Novel Zirconia Phases

- 1) Methyl benzoate
- 2) Ethyl benzoate
- 3) Propyl benzoate
- 4) Butyl benzoate
- 5) Pentyl benzoate



**Column phase:** GROM-SIL Zirconia RP C, 3 µm  
**Column size:** 50 x 2 mm  
**Eluent A:** H<sub>2</sub>O  
**B:** ACN  
**Gradient:** 30 - 80% B (0 - 4.0 min)  
**Flow rate:** 1.0 ml/min  
**Pressure:** 32 MPa  
**Temperature:** RT  
**Detection (UV):** 254 nm (1.2 µl flow cell)  
**Injection:** 2 µl (50-150 µg/ml of each)

GZ RP P ....C0503	Zircon RP P, capillary column,	50 x 0.3 mm
GZ RP P ....C1003	Zircon RP P, capillary column,	100 x 0.3 mm
GZ RP P ....C1503	Zircon RP P, capillary column,	150 x 0.3 mm
GZ RP P ....S0501	Zircon RP P, microbore column,	50 x 1 mm
GZ RP P ....S1001	Zircon RP P, microbore column,	100 x 1 mm
GZ RP P ....S1501	Zircon RP P, microbore column,	150 x 1 mm
GZ RP P ....S0502	Zircon RP P, narrowbore column,	50 x 2 mm
GZ RP P ....S1002	Zircon RP P, narrowbore column,	100 x 2 mm
GZ RP P ....S1502	Zircon RP P, narrowbore column,	150 x 2 mm
GZ RP P ....S0504	Zircon RP P, analytical column,	50 x 4.0 mm
GZ RP P ....S1004	Zircon RP P, analytical column,	100 x 4.0 mm
GZ RP P ....S1504	Zircon RP P, analytical column,	150 x 4.0 mm
GZ RP P ....S0505	Zircon RP P, analytical column,	50 x 4.6 mm
GZ RP P ....S1005	Zircon RP P, analytical column,	100 x 4.6 mm
GZ RP P ....S1505	Zircon RP P, analytical column,	150 x 4.6 mm